

# Widener

## Widener organizer

*instruments not included* 30545



Tray to effectively store the wideners and the handling instruments. The serigraphy allows to immediately identify the widener and the instruments to be used.




## Widener

Nr.		∅ pilot hole	∅	∅	∅	∅	∅	∅
		∅ apical	L. 8	L. 9	L. 10	L. 11	L. 12	L. 13-15
1	15702	2.0	2.20	2.25	2.30	2.35	2.45	2.50
2	15703	2.5	2.60	2.65	2.70	2.75	2.85	2.90
3	15704	2.8	2.95	3.00	3.10	3.15	3.25	3.30
4	15705	3.0	3.15	3.20	3.30	3.35	3.45	3.50
5	15706	3.5	3.60	3.70	3.75	3.85	3.90	4.00
6	15707	4.0	4.15	4.20	4.30	4.35	4.45	4.50



They allow you to gradually enlarge the crest, expanding the available bone and reducing surgical trauma. They increase the transverse volume in presence of thin edentulous ridges with suitable height.



## Insert extension

	21126
<i>o-ring</i>	<i>pack 3pcs</i> 21144



It allows you to use the wideners in less accessible areas or between two dental elements.

